

Laboratory Single Drench Hose

Specification Data Sheet



Accessories

WS-SETAG - Weekly Test Tag For Emergency Shower and/or Eye Wash - 52 Week, Waterproof (Pk 10)

WS-SGN02 - Safety Sign for Emergency Shower - 300 x 225

WS-SGN04 - Safety Sign for Emergency Shower and Eye Wash - 450 x 200

WS-SGN06 - Safety Sign for Emergency Shower and Eye Wash - 300 x 225

To Specify

The SE-920 Drench Hose Emergency Station is designed to provide remarkable flexibility when crisis occurs. Ideal for laboratories and industrial settings, this drench hose can be wall mounted in desired locations and features aerated sprays for soothing relief. A 6 foot, reinforced flexible hose delivers spot on washing to all contaminated areas.

Code:

WS-SE920 - Wall Mounted Laboratory Aerated Eye Wash Hand Held Single Drench Hose, Single Stream.




A family
owned
Australian
business



Performance Standards & Approvals

- Complies to AS4775 - 2007
- All units meet existing ANSI Z358.1 1998 standards and OSHA rules for drench hoses

Typical Illustration



Installation

After unpacking carton, ensure all components are present. Connect water supply to the unit. Turn water supply on and assure all connections are water tight. Activate the unit and allow to flow for 30 seconds to flush line of any impurities.

Warranty

Galvin Engineering products are covered under our Manufacturer's Warranty at www.galvinengineering.com.au. Galvin Engineering products must be installed in accordance with the installation instructions and in accordance with AS3500 and the PCA, relevant Australian Standards and local authorities applicable to product being installed. Water and electrical supply conditions must also comply to the applicable national and/or state standards, failing to comply with these provisions may void the product warranty and affect performance of the product.

Whilst all reasonable care has been taken in compiling the information in this document, the designs, dimensions and products shown are indicative only and should not be relied upon without our prior approval. Due to our policy of continuous development, we reserve the right to alter any details of specifications or products without notice. All information, drawings, technical specifications and product designs remains the property of Galvin Engineering Pty Ltd.

Technical Specifications

- Spray Head: Aerated, ABS Plastic Spray Outlet
- Valve: 10mm Chrome Plated Brass. Self-Closing.
- Activator: Squeeze Handle.
- Flow Pressure: 136 kPa
- Flow Rate (LPM): 6 lpm

Operation

The drench hose is activated by applying pressure to the valve handle. The self-closing valve will open and an aerated flow of water will discharge from the outlet. Water will continue to flow as long as the handle is depressed. Once pressure on the handle is released, the valve will immediately shut off.

Maintenance & Cleaning

Galvin Engineering products should be cleaned with a soft damp cloth using only mild liquid detergent or soap and water. Do not use cleaning agents containing a corrosive acid, scouring agent or solvents. Do not use cream cleaners, as they are abrasive. Use of unsuitable cleaning agents may damage the surface. Any damage caused in this way will not be covered by warranty.

Regular maintenance should be carried out in accordance with the installation and troubleshooting guides provided with the product, relevant Australian Standards, local authorities and state codes.

As part of maintaining a safe and reliable plumbing system, we recommend tapware and valves are operated on at least a daily basis.

